

New & different non-pharmacological mechanism of action

Vibrant is indicated for the treatment of adults with chronic idiopathic constipation who have not experienced relief of their bowel symptoms by using laxative therapies at the recommended dosage for at least one month.



Vibrant developed a new non-pharmacological mechanism of action based on intraluminal mechanical stimulation



1. Activate

Place the capsule in the pod for an automatic activation



2. Ingest

Swallow the capsule with a glass of water



3. Constipation Relief


Intraluminal mechanical stimulation improves natural motility

Efficacy and Safety of the Non-Pharmaceutical Vibrant® System for Chronic Idiopathic Constipation (CIC):

Randomized, double-blind, multicenter, placebo-controlled, phase III trial, Rao SSC et al DDW 2022; Gastroenterology Volume 162, Issue 7, Supplement, S-124

Significantly more CSBMs

Primary Endpoints



Improvement in straining, stool consistency and quality of life

Secondary Endpoints


	Placebo	VC (Vibrant)	P-value
Straining Score	-1.10	-1.66	0.0139
Stool Consistency	0.47	1.00	<0.001
QOL	-0.93	-1.31	0.0020

Safe and well tolerated

	Placebo	VC (Vibrant)
Gastrointestinal disorders- total	9.40%	6.13%
Abdominal pain	3.36%	0.61%
Vomiting	0.67%	1.23%
Nausea	0.67%	1.84%
Abdominal distension	1.34%	0.61%
Diarrhea	-	1.23%

Potential adverse events associated with the use of this device may include abdominal pain, blood in stool and diarrhea. In rare cases, obstruction may occur.

For the full list of side effects, risks and warnings, please scan the QR




or visit: www.vibrantgastro.com/information-for-use/



Get control with our app!

While the treatment works without the app, it makes the experience even better 😊

- E-Diary assists in taking control of the treatment progress
- Get the most out of the treatment with capsule intake reminders
- Easy access to Vibrant support, content & tips



Got more questions?

Scan the QR code for FAQs
Call 866-359-7532 (Mon-Fri 8am-4pm EST),
or visit us at www.vibrantgastro.com



VIBRANT